# UltraBasic

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# Contents

---

Warranty Instructions	4
Disclaimers	5
Meter Specifications	6
Accessories	6
UltraBasic Meter Quick Reference	7
Getting Started	8
Installing and Maintaining Electrodes	10
Standardizing for pH Measurement	11
Using Setup in pH Mode	13
Standardizing for Millivolt Measurement/Relative Millivolts (UP-10 Model Only)	14
Using Setup in mV Mode (UP-10 Model Only)	15
Troubleshooting	16
CE Marking	17

# Warranty Instructions

---

This product is warranted for two years from the date of purchase against defects of material and/or workmanship.

This warranty does not cover batteries and is void if the product is damaged by accident or through unreasonable use, neglect, improper service, or other causes not arising out of defects in material or workmanship.

During the two-year warranty period, if the balance fails to work properly when used under normal conditions we will repair or, at your option, replace the balance free of charge. You must cover the cost of shipping when sending the balance to us; we will bear the shipping costs when returning the balance to you.

Return defective equipment to authorized distributor, or to:

*North and South America:*

**Denver Instrument Company  
6542 Fig Street  
Arvada, Colorado 80004  
USA  
1-800-321-1135  
Tel: 303-431-7255  
Fax: 303-423-4831**

*U.K. and Ireland:*

**Denver Instrument Company  
Denver House, Sovereign Way  
Trafalgar Business Park  
Downham Market  
Norfolk PE38 9SW England  
Tel: 44 136 63862 42  
Fax: 44 136 63862 04**

*Europe, Asia and Australia:*

**Denver Instrument GmbH  
Robert-Bosch-Breite 10  
37079 Goettingen, Germany  
Tel: 49 551 20977 31  
Fax: 49 551 20977 39**

# Disclaimers

---

**Important:**

The operator shall be responsible for any modifications to Denver Instrument Company equipment and for any connections of cables or equipment not supplied by Denver Instrument and must check and if necessary, correct these modifications and connections.

You have purchased a quality precision meter that requires handling with care.

*Please read this manual carefully before operating your new Denver Instrument meter.*

**Caution:**

Changes or modifications not expressly approved by the manufacturer could void the product's warranty.

# Meter Specifications

---

pH	Range	0.00 to 14.00 pH
	Resolution	0.01 pH
	Accuracy	±0.01 pH
mV (UP-10 only)	Range	-1800.0 to 1800.0 mV
	Resolution	0.1 mV
	Accuracy	±0.2 mV or 0.05%, whichever is greater
Temperature	Range	0.0 to 105.0°C
	Resolution	0.1°C
	Accuracy	±0.2°C
Standardization	1, 2 or 3 buffers	
Automatic buffer detection	16 buffers	
	2, 4, 7, 10, 12	
	1, 3, 6, 8, 10, 13	
	1.68, 4.01, 6.86, 9.18, 12.46	

IP-67 rated (waterproof when immersed to 1 meter)  
Auto pH Temperature Compensation with ATC probe  
Automatic pH electrode slope correction for 90-105%

## Accessories

You can order the following accessories for your pH meter:



Replacement UltraBasic pH electrode	<b>301423.1</b>
Plastic-body pH/ATC electrode	<b>300728.1</b>
High-performance glass-body pH/ATC electrode (Tris-compatible)	<b>300729.1</b>
ATC temperature probe	<b>300733.1</b>
Free-standing electrode arm with base	<b>300401.1</b>
Transformer, 115V NA	<b>102077.1</b>

Other electrodes, such as ion-selective electrodes and redox electrodes, are also available. Contact your sales representative for details.

# UltraBasic Meter Quick Reference


---

1. Install a 9-volt battery in the battery compartment or connect an AC adapter to the power connection (optional).
2. Connect an electrode to the Input and ATC connectors (optional).
3. Press **mode** until the desired measurement mode (pH or mV/relative mV) is displayed.
4. Standardize the meter using up to three buffers by immersing the electrode in a buffer, stirring, then pressing **standardize** to enter each buffer.
5. The display shows the current reading in pH, mV, or relative mV units.
6. Press **setup** to review electrode calibration and to clear or select buffer sets.


# Getting Started

## Front Panel Controls


## Digital Display


**Note:** Not all display symbols are used.

## Rear Panel Connectors

## Battery Installation

# Installing and Maintaining Electrodes

---





1. Remove the protective end cover from electrode.
2. Before using your pH electrode, or whenever the electrode is dry, soak overnight in an electrode filling solution, KCl solution or electrode storage solution.
3. Remove the shorting cap on the BNC connector. Install the electrode by plugging it into the **input** connection (push on and turn to lock). An accessory ATC probe can be installed in the **ATC** port for the most accurate measurement.
4. Rinse and blot-dry electrodes between each measurement (**do not wipe**). Rinse electrodes with distilled water or deionized water, or use the next solution to be measured.
5. Store pH electrodes in KCl solution or electrode storage solution. Always leave the filling hole open when in use and closed when in storage. Refill with filling solution when the internal solution level drops to 1" below fill hole.

**Note:** The electrode that comes with the UltraBasic kit is gel-filled. It has no filling solution or filling hole.

# Standardizing for pH Measurement


Because electrodes vary in their response, it is important to standardize your pH meter and electrodes to compensate for electrode variation. The more frequently you standardize, the more accurate your measurements. Standardize **at least** once a day for best results.

**Note:** To exit the standardizing mode, press the “standardize” key again.


1. Immerse electrode in a buffer solution. Stir gently. Allow the electrode reading to stabilize.
2. Press the **mode** key repeatedly until the display indicates the pH mode.
3. Clear existing buffers when performing initial standardization; e.g., for a new set of measurements. Use the **setup** and **enter** keys to clear existing buffers.
4. Press **standardize**. The current buffer set flashes on the display and is detected by the meter. When the signal is stable, or when you press **enter**, the buffer is stored.
5. The meter displays the percent slope of the electrode as 100.0% of the first buffer. On entering a second or third buffer, the meter performs a diagnostic check on the electrode and displays the slope.
6. To enter a second buffer, place the electrode in the second buffer solution, stir, allow time for the electrode to stabilize, and press **standardize** again. The meter detects the buffer.

7. Next, the meter performs a diagnostic test of the electrode. The display indicates the electrode status and the % slope of the electrode.


8. **Error** shown with the 'Slope' symbol indicates that your electrode is not working properly. The electrode response must be between 90 and 105%. Measurements that results in 'Slope Error' are not accepted by the meter. Press **enter** to continue.

9. To enter a third standard, place the electrode in the third buffer solution, stir, allow to stabilize, and press **standardize**. The results are the same as in steps 7 and 8.

10. After all buffers have been entered, the *Standardizing* symbol goes blank and the *Measuring* or *Stable* symbol is displayed to indicate that the meter has returned to *Measuring* mode.

11. Standardize your meter and electrode using at least two buffers with pH values above and below the expected pH of your samples. Stirring with a magnetic stir bar and stirrer provides faster electrode response.

**Note:** If an ATC probe is used, the meter continually adjusts for temperature. In this case, the buffers might vary slightly from the nominal values due to temperature effects.

**Default temperature: 25°C**



**Note:** The temperature is displayed only when an ATC probe is connected.

# Using Setup in pH Mode

---



The **setup** key lets you clear all buffers/standards that have been entered, review standardizing information, and select the desired buffer or set of standards.

**Note:** Press **setup** repeatedly to return to measurement mode.

1. Press **setup**; the word *Clear* flashes on the display. *Continue with this step only if you wish to clear all buffers entered.* To clear all existing buffers, press **enter**. The meter clears all buffers and returns to measuring mode.
2. Press **setup** again to show electrode performance. If two buffers are stored, *OK* is displayed with the electrode symbol, the slope between the first and second buffers, and the two buffer symbols.
3. Press **setup** again to show the electrode slope between the second and third buffers (if three buffers have been entered) and the second and third buffer symbols.



**Note:** The symbols are displayed as described only if buffers are stored.

4. Press **setup** again; *Set* flashes and the first buffer set (U.S.: buffer series) is shown.
5. Press **enter** to select the set of buffers displayed or  
Press **setup** again to show the next set of buffers (Merck buffer set), and again for the third buffer set (NIST buffer set).

## Standardizing for Millivolt Measurement/ Relative Millivolts (UP-10 Model Only)

---


Millivolt measurement is used for electrode diagnostics, titration and measuring redox potential (also called ORP, oxidation reduction potential). ORP measurements indicate the oxidizing or reducing capability of a solution. You can use ORP values to monitor and control solutions requiring a specified amount of oxidants or reductants.


1. Immerse the electrode in a standard solution.
2. Press the **mode** key until your digital display indicates mV mode.
3. Press **standardize** to enter a mV standard and read relative mV.
4. When the signal stabilizes, or when you press **enter**, the current absolute mV value is set to 'zero' (relative millivolts).

## Using Setup in mV Mode (UP-10 Model Only)





---

1. To clear the mV offset and return to absolute millivolt mode, press **setup**. The word *Clear* flashes on the display, and the current relative millivolt offset is shown.
2. To clear the mV standard, press **enter**. The meter returns to absolute mV mode.

# Troubleshooting

1. When the signal from the electrode is out of range, the display shows " - - - ". This may indicate that the electrode is not in solution.
2. When the meter detects an error in the temperature probe, the display shows - - -. If you do not use a temperature probe, the meter applies the default temperature value, 25°C.
3. When an error in electrode response is detected, the display shows *Electrode Error*. During standardization, this message indicates that the electrode exhibits less than 90% or more than 105% of the correct response. The *Electrode Error* message can indicate a bad electrode or bad buffer(s).
4. To test the pH electrode, place it in a good pH 7 buffer. Press **mode** to activate the mV mode, and note the millivolt reading. Make sure the meter is in mV mode and not **relative** mV mode. Repeat the procedure with either a pH 4 or pH 10 buffer. The electrode signal must be within the limits shown below (with the temperature near 25°C).
5. To test the meter for correct operation, install the BNC (input) shorting cap. Press **mode** to select the mV mode, and note the mV reading. Make sure the meter is in mV mode and not relative mV mode. If the meter reads  $0 \pm 0.3$  mV, it is measuring correctly. Note that a long-term drift of 0.1 mV/month since the most recent standardization is specified.

## Electrode Test

pH 7	$0 \pm 30$ mV
pH 4	159 to 186 mV more than pH 7
pH 10	159 to 186 mV less than pH 7

# CE Marking

---

The CE marking affixed to the equipment indicates that the equipment meets the requirements of the following Directive(s) issued by the Council of the European Union:

**Council Directive 89/336/EEC  
“Electromagnetic compatibility  
(EMC)”**

1. Electromagnetic Compatibility
- 1.1 Reference to 89/336/EEC:

Official Journal of the European Communities, No. 2001/C 105/03

**EN 61326-1** Electrical equipment for measurement, control and laboratory use EMC requirements

**Part 1:** General requirements  
Defined immunity to interference: Industrial areas, continuous non-monitored operation  
Limitation of emissions:  
Residential areas, Class B

**Important Note:**

The operator shall be responsible for any modifications to Denver Instrument equipment and for any connections of cables or equipment not supplied by Denver Instrument and must check and, if necessary, correct these modifications and connections. On request, Denver Instrument will provide information on the minimum operating specifications (in accordance with the Standards listed above for defined immunity to interference).

**73/23/EEC “Electrical equipment designed for use within certain voltage limits”**

Applicable European Standards:

**EN 60950** Safety of information technology equipment including electrical business equipment

**EN 61010** Safety requirements for electrical equipment for measurement, control and laboratory use  
Part 1: General requirements

If you use electrical equipment in installations and under ambient conditions requiring higher safety standards, you must comply with the provisions as specified in the applicable regulations for installation in your country.

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Copyright by Denver Instrument GmbH, Goettingen, Germany.  
All rights reserved. No part of this publication may be reprinted  
or translated in any form or by any means without the prior  
written permission of Denver Instrument GmbH.

The status of the information, specifications and illustrations in  
this manual is indicated by the date given below. Denver Instrument  
GmbH reserves the right to make changes to the technology, fea-  
tures, specifications, and design of the equipment without notice.  
Status: May 2003, Denver Instrument GmbH, Goettingen, Germany.

*Europe, Asia and Australia:*

**Denver Instrument GmbH  
Robert-Bosch-Breite 10  
37079 Goettingen, Germany  
Tel: +49-551-20977-30  
Fax: +49-551-20977-39**

*Internet:*

**[www.denverinstrument.com](http://www.denverinstrument.com)**

*U.K. and Ireland:*

**Denver Instrument Company  
Denver House, Sovereign Way  
Trafalgar Business Park  
Downham Market  
Norfolk PE38 9SW England  
Tel: +44-136-63862-42  
Fax: +44-136-63862-04**

*North and South America:*

**Denver Instrument Company  
6542 Fig Street  
Arvada, Colorado 80004  
1-800-321-1135  
Tel: +1-303-431-7255  
Fax: +1-303-423-4831**